**N-Histofine<sup>®</sup>**

**DAB-2V**

Peroxidase Chromogen/Substrate solution

**NICHIREI BIOSCIENCES INC.**

6-19-20, Tsukiji, Chuo-ku, Tokyo 104-8402, JAPAN

Phone: 81-3-3248-2208 Facsimile: 81-3-3248-2243

### 3,3'-DIAMINOBENZIDINE TETRAHYDROCHLORIDE (DAB) SUBSTRATE SOLUTION SUBSTRATE FOR HORSERADISH PEROXIDASE

Package/Code 500 tests/425312F, 1500 tests/425314F

#### I-1. DESCRIPTION

**N-Histofine<sup>®</sup> DAB-2V** is designed as one of reagents to visualize antigen binding primary antibody on tissue sections.

#### I-2. INTENDED USE

**N-Histofine<sup>®</sup> DAB-2V** is designed to allow immunohistochemical tests, to reveal antigens that reacts with a user-supplied primary on tissues and cells. The product is characterized by a coloring reagent/substrate for coloring a detection system reagent labeled with peroxidase. The product is one of the specific chromogen/substrate solutions for peroxidase and develops brown colored precipitate by reaction with peroxidase. The result provides a physiological or pathological state.

For In-vitro Diagnostic Use. For manual staining. Staining results are qualitative. For formalin-fixed paraffin-embedded human tissue sections.

For professional users.

Interpretation must be made within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

#### I-3. PRINCIPLE

The antigen / antibody / Universal Immuno-peroxidase Polymer complex can be prepared by allowing the reagent to react with a primary antibody bound to the antigen on tissue section. The enzymatic activity of this complex results in a colored deposit, thus staining the antigen site.

## II. COMPONENTS

### 1. Reagents provided

Reagent A: DAB solution concentrate (Liquid)	...	[500 tests] 1.2 ml×2	[1500 tests] 1.2 ml×6
3,3'-diaminobenzidine tetrahydrochloride (DAB) chromogen solution	...		
Reagent B: DAB buffer (Liquid)	...	30 ml×2	30 ml×6
Hydrogen peroxide	...		

### 2. Accessories

Tube	...	×1	×1
Chip (For Reagent B, colored)	...	×1	×1
Dropping bottle	...	×1	×1

## III. USAGE

### 1. Preparation of DAB solution

Add 1 drop (approx. 20 µl) of Reagent A Chromogen reagent per 500 µl of Reagent B Chromogen substrate into graduated tube, and place its cap then mix carefully without foaming.

Remove full dosage of the DAB solution from the tube to dropping bottle using pipette as well for use.

Store in light shielded condition at 2 to 8°C and use within 2 weeks after the preparation.

- Note:** 1. Use DAB solution removing foams or exclusive of foams when its surface is foamed by mixing.  
2. Do not use mixture of DAB solution prepared different day.

### 2. Addition and Reaction of DAB solution

1. Apply DAB solution using dropping bottle covering the tissue section on the slide completely and incubate at room temperature for 5 to 20 minutes after addition, reaction and washing of peroxidase conjugate in the immuno- histochemical staining process.
2. After adequate chromogenic reaction, wash the slide well with tap water.
3. Mount the slide with permanent mounting media and fix with cover slip after dehydrated in graded series of alcohol and cleared in xylene.

**Note:** 1. Adjust the reaction time by microscopic observation.

2. At the same time, staining of negative control slide using such as PBS in place of primary antibody is recommended to be prepared for checking the staining result.

### 3. Accessories

#### Tube

For mixing of Reagent A and Reagent B. (Refer to aforementioned Preparation of DAB solution.)

#### Chip (For Reagent B, colored)

Reagent B is provided with translucent chip placed on its bottle. However the chip is not allowed for applying solution drop by drop, therefore it is not suitable for small volume preparation. For such preparation, replace the translucent chip by colored chip provided exclusively for dropping is recommended. The colored chip is suitable for applying solution drop by drop.

#### Dropping bottle

Dropping bottle is used for storage of reagents, not only for dropping of DAB solution.

**Note:** Tube and dropping bottle are reusable after wash with purified water and dry completely.

### IV. STAINING PROCEDURE

Refer to the instruction\* of **M**-Histofine® Simple Stain MAX PO (MULTI) for staining procedures such as Principle, Materials required but not supplied, Specimen preparation, Control slides, Interpretation of results and General limitation. Validation is required when other detection systems besides above are used.

\*Download of the instruction is available at following website:

[http://www.nichirei.co.jp/bio/english/products/pdf/CE\\_MAXPO\\_MULTI.pdf](http://www.nichirei.co.jp/bio/english/products/pdf/CE_MAXPO_MULTI.pdf)

### V. STORAGE, HANDLING & SHELF LIFE

1. Store at 2-8°C.
2. Store prepared reagent, DAB solution, in light shielded condition at 2 to 8°C and use within 2 weeks after the preparation.
3. Confirm expiration date indicated on the labels of vial and carton before use.
4. Adjust this product to room temperature before use.
5. Refrigerate this product promptly after use.
6. The reagent is stable 18 months after manufacturing.

### VI. PRECAUTIONS

1. Before using this reagent, read these instructions.
2. Do not use reagents after the expiration date.
3. For professional users.
4. Specimens, before and after fixation, and all other materials exposed to them, should be handled like biohazardous materials with proper precautions.
5. Inhalation or ingestion of the highly allergic fixative formaldehyde is harmful. Wear protective mask. If swallowed, induce vomiting. If skin or eye contact occurs, wash thoroughly with water.
6. Organic reagents are flammable. Do not use near open flame.
7. Never pipette reagents by mouth and avoid their contact with skin, mucous membranes and clothes.
8. Avoid microbial contamination of reagents as incorrect result may occur.
9. Avoid splashing of reagents or generation of aerosols.
10. As a chromogen, DAB solution should be handled with care for it contains carcinogen.
11. Unused solution should be disposed of according to local, State and Federal regulations.
12. Any serious incident that has occurred in relation to **M**-Histofine® DAB-2V shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

### VII. REFERENCE

- (1) Mokry, et al. Versatility of immunohistochemical reactions: comprehensive survey of detection systems. Acta medica 1996 Apr 39:129-40
- (2) Yamada K, et al. In vitro assessment of antitumor immune responses using tumor antigen proteins produced by transgenic silkworms. J Mater Sci Mater Med. 2021 May 17;32(6):58.

 <b>REF</b>	Catalog/Code Number		Temperature Limitations	 <b>IVD</b>	In Vitro Diagnostic Medical Device
	Manufacturer	 <b>LOT</b>	Batch Code		Contains Sufficient for <N> Tests
	Use By		Consult Instructions for Use	 <b>EC REP</b>	Authorized Representative in the European Community
 <b>CE</b>	CE-mark, code of the notified body		For IVD Performance Evaluation only	 <b>SAMPLE</b>	Sample

**NICHIREI BIOSCIENCES INC.**  
6-19-20, Tsukiji, Chuo-ku, Tokyo 104-8402, JAPAN  
Phone: 81-3-3248-2208, Facsimile: 81-3-3248-2243

MedEnvoy Global B.V.  
Prinses Margrietplantsoen 33-Suite 123  
2595 AM The Hague  
The Netherlands

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### V. CONSERVATION, MANIPULATION & PEREMPTION

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 REF	Catalog/Code Number		Temperature Limitations	 IVD	In Vitro Diagnostic Medical Device
	Manufacturer	 LOT	Batch Code		Contains Sufficient for (N) Tests
	Use By		Consult Instructions for Use	 EC REP	Authorized Representative in the European Community
	CE-mark, code of the notified body		For IVD Performance Evaluation only	 SAMPLE	Sample

NICHIREI BIOSCIENCES INC.  
6-19-20, Tsukiji, Chuo-ku, Tokyo 104-8402, JAPAN  
Phone: 81-3-3248-2208, Facsimile: 81-3-3248-2243

MedEnvoy Global B.V.  
Prinses Margrietplantsoen 33-Suite 123  
2595 AM The Hague  
The Netherlands

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