



2 Step Polymer Detection System


N-Histofine[®] High Stain[™] HRP (MULTI)

Feature

N-Histofine[®] High Stain[™] HRP (MULTI) is a **2** Step Polymer Detection System to IHC staining providing more amplified staining intensity compared with conventional 1 step polymer detection system. This system is applicable to both of mouse and rabbit primary antibodies and is for formalin-fixed, paraffin-embedded tissue sections.

The Universal Immuno-enzyme Polymer (UIP) method* developed by Nichirei Biosciences Inc. was adopted for this system.

* US. Patent No.6252053, 6613564 and EPC Patent No.0992794.



NICHIREI BIOSCIENCES INC.




- High intensity of staining
- Low expression of antigen is detectable
- No background unaffected by endogenous biotin
- Applicable to mouse or rabbit primary antibodies
- Shortened reaction time only for 20 minutes




Saving your primary antibodies

CDX-2 Rabbit-Monoclonal antibody (Code: 418011F)


Human colon cancer Positive reaction is observed in nuclei of sporadic tumor cells.



Competitive Product



Cyclin D1 Rabbit-Monoclonal antibody (Clone: SP4) (Code: 413521F)



EGFR Mouse-Monoclonal antibody (Clone: 31G7) (Code: 423701F)

Human colon cancer Positive reaction is observed in cytomembrane and cytoplasm of sporadic tumor cells.

Human breast cancer (metastasis to liver)


 M-Histofine® High Stain™ HRP (MULTI)
 Competitive Product

 Dilution of primary antibody
 Image: Competitive Product

 × 1
 Image: Competitive Product

 × 5
 Image: Competitive Product

Principle & Procedure



Product List

Liquid. Ready to us				
Product Name	Tests	Volume	Code	For use with
M-Histofine [®] High Stain [™] HRP (MULTI) High Stain [™] Solution A High Stain [™] Solution B	170	17ml each	414481F	Mouse and Rabbit primary antibody
	1,000	17ml x 6 each	414483F	

NICHIREI BIOSCIENCES INC.

6-19-20, TSUKIJI, CHUO-KU, TOKYO, 104-8402, JAPAN PHONE: 81-3-3248-2207 FAX: 81-3-3248-2243 http://www.nichirei.co.jp/bio/english