

# IHC staining system

## For detection of lower expression of ALK fusion proteins



# **N-Histofine<sup>®</sup> ALK Detection KIT** Code: 417071F Size: 20 tests

**M**-Histofine<sup>®</sup> ALK Detection Kit detects anaplastic lymphoma kinase (ALK) proteins in tumor cells in paraffin-embedded tissue specimens by IHC staining and determines presence of such protein expression.

### Background

The *ALK* gene was identified in 1994 as a gene fused to the nucleophosmin (*NPM*) gene in anaplastic large-cell lymphoma (ALCL) with t(2;5)(p23;q35) translocation. This gene is located at 2p23 and encodes for a receptor-type tyrosine kinase, which belongs to the insulin receptor family. The ALK protein have a kinase domain in its intracellular domain, and its function is associated with the promotion of cell growth and inhibition of apoptosis.


Subsequently, the ALK gene has been reported to form ALK fusion genes fused with ATIC, CLTC, MSN, TPM3, TPM4, TFG, MYH9 and ALO17 genes in ALCL and also ATIC, CARS, CLTC, DCTN1, TPM3, TPM4, PPFIBP1, RANBP2 and SEC31L1 genes in inflammatory myofibroblastic tumor (IMT).


The proteins produced from these fusion genes are constantly activated by forming dimers and led to cancerous change.


Recently, other *ALK* fusion genes with *EML4* gene, *KIF5B* gene or *KCL1* gene in non-small cell carcinoma of lung, *SEC31A* gene or *SQSTM1* gene in ALK-positive large B-cell lymphoma and *VCL* gene in renal cell cancer have been also reported.

### **Staining Image**

Strong Positive Image







Lung adenocarcinoma: Positive for cytoplasm of tumor cells. Weak to strong positive of staining levels regarding the expression level of ALK fusion protein. Hot bath treatment (+)

## Kit Components

	Product Name		Size	Code
<b>N</b> -	Histofine <sup>®</sup> ALK Detection	n Kit	20 tests	417071F
Vial No.	Components	Constitue	ents	Volume
1	Blocking Reagent	3 V/V% Hydrogen	peroxide	4ml × 1
2	Primary Antibody anti-ALK mouse monoclonal antibody (5A-		A4) 2ml × 1	
3	Negative Control Mouse IgG		2ml × 1	
4	Bridge Reagent			4ml × 1
5	Peroxidase Labeled Empower Reagen	nt		4ml × 1
6	Chromogen Substrate	3,3'-Diaminobenzid	line tetrahydrochloride	e 0.5ml × 1
7	Substrate Buffer Solution			0.5ml × 1
8	Chromogen Reagent	0.6 V/V% Hydroger	0.6 V/V% Hydrogen peroxide solution	
9	ALK Antigen Retrieval Solution A 150ml × 1			150ml × 1
10	ALK Antigen Retrieval Solution B			150ml × 1
Product Name			Size	Code
<b>M</b> -Histofine <sup>®</sup> ALK Control Slides			5 slides	417081F

#### Notes for Determination

Due to this kit detects ALK proteins, ALK fusion proteins as well as full-length ALK protein are reacted. Consequently, results of slight positive to positive for tumors\* that infrequently express full-length ALK protein are observed. However, discrimination between ALK fusion proteins and full-length ALK protein is not available. Therefore, considering expression possibility of ALK fusion proteins, confirmation of the presence or absence of *ALK* fusion genes by using FISH method is preferable in this regard.

- \* Large-cell neuroendocrine carcinomas of the lung, small-cell lung carcinomas and rhabdomyosarcomas (particularly alveolar rhabdomyosarcomas).
  - Weak positive Images in large cell neuroendocrine carcinoma of the lung

Weak positive staining is observed in the specimen which expresses full-length ALK protein. Hot bath treatment (+)

### References


- (1) Shiota M., et al: Hyperphosphorylation of a novel 80 kDa protein-tyrosine kinase similar to Ltk in a human Ki-1 lymphoma cell line, AMS3. Oncogene 9: 1567-1574, 1994
- (2) Morris SW., et al: Fusion of a kinase gene, ALK, to a nucleolar protein gene, NPM, in non-Hodgkin's lymphoma. Science. 263: 1281-1284, 1994
- (3) Takeuchi K., et al: KIF5B-ALK, a novel fusion oncokinase identified by an immunohistochemistry-based diagnostic system for ALK-positive lung cancer. Clin Cancer Res. 15(9): 3143-3149, 2009
- (4) Takeuchi K., et al: Pulmonary inflammatory myofibroblastic tumor expressing a novel fusion, PPFIBP1-ALK: reappraisal of anti-ALK immunohistochemistry as a tool for novel ALK-fusion identification. Clin Cancer Res. 17(10): 3341-3348, 2011
- (5) Sugawara E., et al. Identification of Anaplastic Lymphoma Kinase Fusions in Renal Cancer. Cancer, 2012


Manufacturer

## NICHIREI BIOSCIENCES INC.

6-19-20, TSUKIJI, CHUO-KU, TOKYO 104-8402, JAPAN PHONE : 81-3-3248-2207 FAX : 81-3-3248-2243 http://www.nichirei.co.jp/bio/english

## **N**-Histofine<sup>®</sup> ALK Detection KIT






## N-Histofine<sup>®</sup> ALK Control Slides Code: 417081F Size: 5 slides


**N**-Histofine® ALK Control Slides are used as a standard for discrimination of positive tissue blocks possess expressed ALK protein.

This product is used for validation of reagent performance and staining technique in the IHC staining (including cytology staining) with **N**-Histofine® ALK Detection Kit.

Both formalin fixed paraffin embedded cell lines of NCI-H2228 (Positive) and SK-BR-3 (Negative) are mounted on each slide.



No stain is observed in the cytoplasm of all cells. Treatment with hot bath (+)



Strong stains are observed in the cytoplasm of the majority of cells. Treatment with hot bath (+)  $% \left( 1-\frac{1}{2}\right) =0$ 

