NICHIREI

N-Histofine[®]
DAB-3S kit

Peroxidase Chromogen/Substrate kit

500 tests Code: 415192F

1500 tests Code: 415194F

3, 3'-DIAMINOBENZIDINE TETRAHYDROCHLORIDE (DAB) SUBSTRATE KIT
SUBSTRATE FOR HORSE RADISH PEROXIDASE

I-1. DESCRIPTION

N-Histofine[®] DAB 3S kit is designed as one of reagents to visualize antigen binding primary antibody on tissue sections.

I-2. INTENDED USE

N-Histofine[®] DAB 3S kit is designed to allow immunohistochemical tests, to reveal antigens that reacts with a user-supplied primary on tissues and cells. The product is characterized by a coloring reagent/substrate for coloring a detection system reagent labeled with peroxidase. The product is one of the specific chromogen/substrate solutions for peroxidase and develops brown colored precipitate by reaction with peroxidase. The result provides a physiological or pathological state.

For In-vitro Diagnostic Use. For manual staining. Staining Results are qualitative. For formalin-fixed paraffin-embedded human tissue sections.

For professional users.

Interpretation must be made within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

I-3. PRINCIPLE

The antigen / antibody / Universal Immuno-peroxidase Polymer complex can be prepared by allowing the reagent to react with a primary antibody bound to the antigen on tissue section. The enzymatic activity of this complex results in a colored deposit, thus staining the antigen site.

II. FORM

Liquid.

This product consists of followings:

Reagent A: DAB solution concentrate	...	[500 tests] 3 ml × 1 bottle	[1500 tests] 9 ml × 1 bottle
Reagent B: Enhancer solution concentrate	...	3 ml × 1 bottle	9 ml × 1 bottle
Reagent C: Hydrogen peroxide solution concentrate	...	3 ml × 1 bottle	9 ml × 1 bottle
Dropping bottle	...	1 bottle	1 bottle
Tube	...	1 tube	1 tube

III. DESCRIPTION

N-Histofine[®] DAB-3S kit is the chromogen/substrate system suitable for use in peroxidase-based immunohistochemical staining. DAB may produce brown precipitates at the site of the target antigen or nucleic acid reacting with peroxidase.

IV. PREPARATION of Chromogen/Substrate solution

1. Add 1 drop of Reagent A (DAB solution concentrate) to 1 ml of distilled water in the tube, and mix well.
2. Add 1 drop of Reagent B (Enhancer solution concentrate) to 1 ml of above diluted Reagent A solution, and mix well.
3. Add 1 drop of Reagent C (Hydrogen peroxide solution concentrate) to 1 ml of above mixture solution of Reagent A and Reagent B, and mix well as final mixture*¹.

*¹: Final mixture is stable for 14 days at 2-8°C under the light shielded condition.

Precipitates may be produced with time. The precipitates adhered on slides or tissues will not affect on the staining, however, in some cases, obtaining clear image may not be suitable.

To avoid such adhesion, following procedures are recommended.

1. Use supernatant solution of the mixture.
2. Rinse the slides in PBS containing 0.1% Tween 20 or TBS containing Tween 20 at room temperature for 3 minutes, at respective steps of before and after the staining described in USAGE.

V. USAGE

Adjust the final mixture to room temperature prior to use.

React the objective tissue section with peroxidase conjugate in advance.

At the same time, staining of Negative control slide using PBS in place of primary antibody is recommended to be prepared for checking the staining result.

For color development, incubate the tissue section on the slide with 1-2 drops of the final mixture using Dropping bottle*2 at room temperature for 5-20 minutes.

Adjust the reaction time by microscopic observation by contrast with the Negative control slide.

After adequate coloration, wash the slide with tap water.

Mount the slide with permanent mounting media and fix with cover slip after washing, dehydrated in graded series of alcohol and cleared in xylene.

- *2: 1) Dropping bottle is used for dropping the final mixture, not for mixing or storage of reagents.
2) When the final mixture is removed from tube to the dropping bottle, adequate pipetting is required.
3) Rinse the dropping bottle thoroughly with distilled water and dry after use.

VI. STORAGE, HANDLING & SHELF LIFE

1. Store at 2-8°C.
2. Confirm expiration date indicated on the labels of vial and carton before use.
3. Refrigerate this product promptly after use.
4. The reagent is stable 18 months after manufacturing.

VII. PRECAUTIONS

1. Before using this reagent, read these instructions.
2. Do not use reagents after the expiration date.
3. For professional users.
4. Specimens, before and after fixation, and all other materials exposed to them, should be handled like biohazardous materials with proper precautions.
5. Inhalation or ingestion of the highly allergic fixative formaldehyde is harmful. Wear protective mask. If swallowed, induce vomiting. If skin or eye contact occurs, wash thoroughly with water.
6. Organic reagents are flammable. Do not use near open flame.
7. Never pipette reagents by mouth and avoid their contact with skin, mucous membranes and clothes.
8. Avoid microbial contamination of reagents as incorrect result may occur.
9. Avoid splashing of reagents or generation of aerosols.
10. As a chromogen, DAB solution should be handled with care for it contains carcinogen.
11. Unused solution should be disposed of according to local, State and Federal regulations.
12. Any serious incident that has occurred in relation to **N**-Histofine® DAB-3S kit shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

VIII. REFERENCE

- (1) Umayahara T, et al. Protective role of Galectin-7 for skin barrier impairment in atopic dermatitis. Clin Exp Allergy. 2020 Aug;50(8):922-931.
- (2) Skálová H, et al. mpact of chemotherapy on the expression of claudins and cadherins in invasive breast cancer. Exp Ther Med. 2019 Oct;18(4):3014-3024.

 REF	Catalog/Code Number		Temperature Limitations	 IVD	In Vitro Diagnostic Medical Device
	Manufacturer	 LOT	Batch Code		Contains Sufficient for <N> Tests
	Use By		Consult Instructions for Use	 EC REP	Authorized Representative in the European Community
	CE-mark, code of the notified body		For IVD Performance Evaluation only	 SAMPLE	Sample

NICHIREI BIOSCIENCES INC.
6-19-20, Tsukiji, Chuo-ku, Tokyo 104-8402, JAPAN
Phone:81-3-3248-2208, Facsimile:81-3-3248-2243

MedEnvoy Global B.V.
Prinses Margrietplantsoen 33-Suite 123
2595 AM The Hague
The Netherlands

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