






# **N**-Histofine<sup>®</sup> DAB-2V

## DAB Substrate kit for HRP Detection

**N**-Histofine<sup>®</sup> DAB-2V is a detection system for immunohistochemical staining to detect HRP.

- Simple step for preparation of substrate solution.
- Higher sensitivity than that of other DAB substrate kits.
- Prepared substrate solution is available for two weeks.




· · · · · · · · · · · · · · · · · · ·			
Reagent A: DAB solution con Reagent B: DAB buffer	centrate 1.2 ml x2 30 ml x2	1.2 ml x6 30 ml x6	CE
Tube Chip (For Reagent B)	x1 x1	x1 x1	
Dropping bottle	x1	x1	

## NICHIREI BIOSCIENCES INC.

### **N**-Histofine<sup>®</sup> DAB-2V

#### Human tonsil cancer


Positive reaction is observed in cytomembrane of sporadic tumor cells.




## Actin, Smooth Muscle Mouse Monoclonal antibody (Code: 412021F)

### **N**-Histofine<sup>®</sup> DAB-2V

Human small intestine Positive reaction is observed in cytoplasm of sporadic smooth muscle cells.



Company A



## NICHIREI BIOSCIENCES INC.

6-19-20, TSUKIJI, CHUO-KU, TOKYO, 104-8402, JAPAN PHONE: 81-3-3248-2207 FAX: 81-3-3248-2243 http://www.nichirei.co.jp/bio/english